



Comprehensive Tumor Profiling for New York State


A better foundation for Molecular Intelligence

The Caris Molecular Intelligence® comprehensive tumor profiling approach to assess DNA, RNA and proteins reveals a molecular blueprint to guide more precise and individualized treatment decisions from among 60+ FDA-approved therapies.



DNA

Mutations & Indels



RNA

Whole Transcriptome Sequencing
(Fusions & Variant Transcripts)



Protein

Immunohistochemistry

Technical Specifications

Sufficient tumor must be present to complete all analysis. If you have any questions, please contact Customer Support at (888) 979-8669.

| Technical Information | IHC | CISH | FISH |
|---|--|---|---|
| Sample Requirements <i>(see requisition for full details)</i> | 1 unstained slide at 4µm thickness from FFPE block, with evaluable tumor present, per IHC test | 1 unstained slide at 4µm thickness from FFPE block, with at least 20-100 evaluable tumor cells present, per CISH test | 2 unstained slides at 4µm thickness from FFPE block, with at least 100 evaluable cells present and 10% tumor, per FISH test |
| Sensitivity/Specificity | >95% | >95% | >95% |

| Technical Information | Next-Generation Sequencing (DNA) | Whole Transcriptome Sequencing (RNA) |
|---|--|--|
| Sample Requirements | FFPE block or 10 unstained slides with a minimum of 20% malignant origin. Needle biopsy is also acceptable (4-6 cores). | FFPE block or 2-5 unstained slides with a minimum of 20% malignant origin. Needle biopsy is also acceptable (4-6 cores). |
| Tumor Enrichment (when necessary) | Microdissection to isolate and increase the number of cancer cells to improve test performance and increase the chance for successful testing from small tumor samples | |
| Number of Genes | 592 genes | ~22,000 genes |
| Average Depth of Coverage (DNA) Average Read Count (RNA) | >750X | 60 million |
| Positive Percent Agreement (PPA) | > 95% for base substitutions at ≥ 5% mutant allele frequency; > 95% for indels at ≥ 5% mutant allele frequency; >90% for copy number alterations (amplifications ≥ 6 copies) | >97% |
| Negative Percent Agreement (NPA) | >99% | >99% |
| Genomic Signatures | Microsatellite Instability (MSI), Tumor Mutational Burden (TMB) | - |

Biomarker Analysis by Tumor Type for New York State

The information below details the biomarkers analyzed by technology for the tumor type submitted. Before ordering testing services, please refer to the profile menu online (www.CarisMolecularIntelligence.com/profiling-menu) to view the most up-to-date listing of biomarkers that will be performed. Tests may vary if insufficient tumor samples are submitted.

| MI Profile™ | | | | | |
|------------------------------------|--|--|--------------------|--------------------------------------|--|
| Tumor Type | Immunohistochemistry (IHC) | Next-Generation Sequencing (NGS) <i>(see reverse for gene list)</i> | | Whole Transcriptome Sequencing (WTS) | Other |
| | | DNA | Genomic Signatures | RNA | |
| Bladder | MMR, PD-L1 (SP142 and 22c3) | Mutation, Indels | MSI | Fusion Analysis | |
| Breast | AR, ER, Her2/Neu, MMR, PD-L1 (SP142), PR, PTEN | Mutation, Indels | MSI | Fusion Analysis | TOP2A (<i>Chromogenic in situ Hybridization</i>) |
| Cancer of Unknown Primary - Female | AR, ER, Her2/Neu, MMR, PD-L1(SP142) | Mutation, Indels | MSI | Fusion Analysis | |
| Cancer of Unknown Primary - Male | AR, Her2/Neu, MMR, PD-L1 (SP142) | Mutation, Indels | MSI | Fusion Analysis | |
| Cervical | ER, MMR, PD-L1 (22c3), PR | Mutation, Indels | MSI | Fusion Analysis | |
| Cholangiocarcinoma/Hepatobiliary | Her2/Neu, MMR, PD-L1 (SP142) | Mutation, Indels | MSI | Fusion Analysis | Her2 (<i>Chromogenic in situ Hybridization</i>) |
| Colorectal and Small Intestinal | Her2/Neu, MMR, PD-L1 (SP142), PTEN | Mutation, Indels | MSI | Fusion Analysis | |
| Endometrial | ER, MMR, PD-L1 (SP142), PR, PTEN | Mutation, Indels | MSI | Fusion Analysis | |
| Esophageal Cancer | Her2/Neu, MMR, PD-L1 (22c3) | Mutation, Indels | MSI | Fusion Analysis | |
| Gastric/GEJ | Her2/Neu, MMR, PD-L1 (22c3) | Mutation, Indels | MSI | Fusion Analysis | Her2 (<i>Chromogenic in situ Hybridization</i>) |
| GIST | MMR, PD-L1 (SP142), PTEN | Mutation, Indels | MSI | Fusion Analysis | |
| Glioma | MMR, PD-L1 (SP142) | Mutation, Indels | MSI | Fusion Analysis | MGMT Methylation (<i>Pyrosequencing</i>) |
| Head & Neck | MMR, p16, PD-L1 (22c3) | Mutation, Indels | MSI | Fusion Analysis | HPV (<i>Chromogenic in situ Hybridization</i>), reflex to confirm p16 result |
| Kidney | MMR, PD-L1 (SP142) | Mutation, Indels | MSI | Fusion Analysis | |
| Melanoma | MMR, PD-L1 (SP142) | Mutation, Indels | MSI | Fusion Analysis | |
| Merkel Cell | MMR, PD-L1 (SP142) | Mutation, Indels | MSI | Fusion Analysis | |
| Neuroendocrine/Small Cell Lung | MMR, PD-L1 (SP142) | Mutation, Indels | MSI | Fusion Analysis | |
| Non-Small Cell Lung | ALK, MMR, PD-L1 (22c3), PTEN | Mutation, Indels | MSI | Fusion Analysis | |
| Ovarian | ER, MMR, PD-L1 (SP142), PR | Mutation, Indels | MSI | Fusion Analysis | |
| Pancreatic | MMR, PD-L1 (SP142) | Mutation, Indels | MSI | Fusion Analysis | |
| Prostate | AR, MMR, PD-L1 (SP142) | Mutation, Indels | MSI | Fusion Analysis | |
| Salivary Gland | AR, Her2/Neu, MMR, PD-L1 (SP142) | Mutation, Indels | MSI | Fusion Analysis | |
| Sarcoma | MMR, PD-L1 (SP142) | Mutation, Indels | MSI | Fusion Analysis | |
| Thyroid | MMR, PD-L1 (SP142) | Mutation, Indels | MSI | Fusion Analysis | |
| Uterine Serous | ER, Her2/Neu, MMR, PD-L1 (SP142), PR, PTEN | Mutation, Indels | MSI | Fusion Analysis | Her2 (<i>Chromogenic in situ Hybridization</i>) |
| Vulvar Cancer (SCC) | ER, MMR, PD-L1 (22c3), PR | Mutation, Indels | MSI | Fusion Analysis | |
| Other Tumors | MMR, PD-L1 (SP142) | Mutation, Indels | MSI | Fusion Analysis | |

MMR = Mismatch Repair proteins: MLH1, MSH2, MSH6, PMS2

For PD-L1 IHC testing, the antibody tested is listed above. For non-urothelial bladder cancers, PD-L1 clones SP142 and 22c3 are performed.

Next-Generation Sequencing Gene List

| Next-Generation Sequencing – Genomic Stability Testing (DNA) | | | | | | | | | |
|---|--------------|------------------|-----------|-----------------|--------------------------------|-----------------|----------|---------------------------|----------|
| Microsatellite Instability (MSI) | | | | | Tumor Mutational Burden (TMB)* | | | | |
| Next-Generation Sequencing – Point Mutations and Indels (DNA) | | | | | | | | | |
| ABI1 | BRD4 | CRLF2 | FOXO4 | HOXC11 | KLF4 | MUC1 | PAK3 | RHOH | TAL2 |
| ABL1 | BTG1 | DDB2 | FSTL3 | HOXC13 | KLK2 | MUTYH | PATZ1 | RNF213 | TBL1XR1 |
| ACKR3 | BTK | DDIT3 | GATA1 | HOXD11 | LASP1 | MYCL (MYCL1) | PAX8 | RPL10 | TCEA1 |
| AKT1 | C15orf65 | DNM2 | GATA2 | HOXD13 | LMO1 | NBN | PDE4DIP | SEPT5 | TCL1A |
| AMER1 (FAM123B) | CBLC | DNMT3A | GNA11 | HRAS | LMO2 | NDRG1 | PHF6 | SEPT6 | TERT |
| AR | CD79B | EIF4A2 | GPC3 | IKBKE | MAFB | NKX2-1 | PHOX2B | SFPQ | TFE3 |
| ARAF | CDH1 | ELF4 | HEY1 | INHBA | MAX | NONO | PIK3CG | SLC45A3 | TFPT |
| ATP2B3 | CDK12 | ELN | HIST1H3B | IRS2 | MECOM | NOTCH1 | PLAG1 | SMARCA4 | THRAP3 |
| ATRX | CDKN2B | ERCC1 | HIST1H4I | JUN | MED12 | NRAS | PMS1 | SOC31 | TLX3 |
| BCL11B | CDKN2C | ETV4 | HLF | KAT6A (MYST3) | MKL1 | NUMA1 | POU5F1 | SOX2 | TMPRSS2 |
| BCL2 | CEBPA | FAM46C | HMG2N2P46 | KAT6B | MLLT11 | NUTM2B | PPP2R1A | SPOP | UBR5 |
| BCL2L2 | CHCHD7 | FANCF | HNF1A | KCNJ5 | MN1 | OLIG2 | PRF1 | SRC | VHL |
| BCOR | CNOT3 | FEV | HOXA11 | KDM5C | MPL | OMD | PRKDC | SSX1 | WAS |
| BCORL1 | COL1A1 | FOXL2 | HOXA13 | KDM6A | MSN | P2RY8 | RAD21 | STAG2 | ZBTB16 |
| BRD3 | COX6C | FOXO3 | HOXA9 | KDSR | MTCP1 | PAFAH1B2 | RECQL4 | TAL1 | ZRSR2 |
| Next-Generation Sequencing – Point Mutations, Indels and Copy Number Alterations* (DNA) | | | | | | | | | |
| ABL2 | BRCA1 | CREB3L1 | ETV1 | GAS7 | KMT2A (MLL) | MYCN | PER1 | RUNX1 | TFEB |
| ACSL3 | BRCA2 | CREB3L2 | ETV5 | GATA3 | KMT2C (MLL3) | MYD88 | PICALM | RUNX1T1 | TFG |
| ACSL6 | BRIP1 | CREBBP | ETV6 | GID4 (C17orf39) | KMT2D (MLL2) | MYH11 | PIK3CA | SBDS | TFRC |
| ADGRA2 | BUB1B | CRKL | EWSR1 | GMP5 | KNL1 | MYH9 | PIK3R1 | SDC4 | TGFBR2 |
| AFDN | CACNA1D | CRTC1 | EXT1 | GNA13 | KRAS | NACA | PIK3R2 | SDHAF2 | TLX1 |
| AFF1 | CALR | CRTC3 | EXT2 | GNAQ | KTN1 | NCKIPSD | PIM1 | SDHB | TNFAIP3 |
| AFF3 | CAMTA1 | CSF1R | EZH2 | GNA5 | LCK | NCOA1 | PML | SDHC | TNFRSF14 |
| AFF4 | CANT1 | CSF3R | EZR | GOLGA5 | LCP1 | NCOA2 | PMS2 | SDHD | TNFRSF17 |
| AKAP9 | CARD11 | CTCF | FANCA | GOPC | LGR5 | NCOA4 | POLE | SEPT9 | TOP1 |
| AKT2 | CARS | CTLA4 | FANCC | GPHN | LHFPL6 | NF1 | POT1 | SET | TP53 |
| AKT3 | CASP8 | CTNNA1 | FANCD2 | GRIN2A | LIFR | NF2 | POU2AF1 | SETBP1 | TPM3 |
| ALDH2 | CBFA2T3 | CTNNA1 | FANCE | GSK3B | LPP | NFE2L2 | PPARG | SETD2 | TPM4 |
| ALK | CBFB | CYLD | FANCG | H3F3A | LRIG3 | NFIB | PRCC | SF3B1 | TPR |
| APC | CBL | CYP2D6 | FANCL | H3F3B | LRP1B | NFKB2 | PRDM1 | SH2B3 | TRAF7 |
| ARFRP1 | CBLB | DAXX | FA5 | HERPUD1 | LYL1 | NFKBIA | PRDM16 | SH3GL1 | TRIM26 |
| ARHGAP26 | CCDC6 | DDR2 | FBXO11 | HGF | MAF | NIN | PRKAR1A | SLC34A2 | TRIM27 |
| ARHGGEF12 | CCNB1IP1 | DDX10 | FBXW7 | HIP1 | MALT1 | NOTCH2 | PRRX1 | SMAD2 | TRIM33 |
| ARID1A | CCND1 | DDX5 | FCRL4 | HMG1 | MAML2 | NPM1 | PSIP1 | SMAD4 | TRIP11 |
| ARID2 | CCND2 | DDX6 | FGF10 | HMG2 | MAP2K1 (MEK1) | NR4A3 | PTCH1 | SMARCB1 | TRRAP |
| ARNT | CCND3 | DEK | FGF14 | HNRNP2B1 | MAP2K2 (MEK2) | NSD1 | PTEN | SMARCE1 | TSC1 |
| ASPSR1 | CCNE1 | DICER1 | FGF19 | HOOX3 | MAP2K4 | NSD2 | PTPN11 | SMO | TSC2 |
| ASXL1 | CD274 (PDL1) | DOT1L | FGF23 | HSP90AA1 | MAP3K1 | NSD3 | PTPRC | SNX29 | TSHR |
| ATF1 | CD74 | EBF1 | FGF3 | HSP90AB1 | MCL1 | NT5C2 | RABEP1 | SOX10 | TTL |
| ATIC | CD79A | ECT2L | FGF4 | IDH1 | MDM2 | NTRK1 | RAC1 | SPECC1 | U2AF1 |
| ATM | CDC73 | EGFR | FGF6 | IDH2 | MDM4 | NTRK2 | RAD50 | SPEN | USP6 |
| ATP1A1 | CDH11 | ELK4 | FGFR1 | IGF1R | MDS2 | NTRK3 | RAD51 | SRGAP3 | VEGFA |
| ATR | CDK4 | ELL | FGFR1OP | IKZF1 | MEF2B | NUP214 | RAD51B | SRSF2 | VEGFB |
| AURKA | CDK6 | EML4 | FGFR2 | IL2 | MEN1 | NUP93 | RAF1 | SRSF3 | VTG1A |
| AURKB | CDK8 | EMSY | FGFR3 | IL21R | MET | NUP98 | RALGDS | SS18 | WDCP |
| AXIN1 | CDKN1B | FGFR4 | FGFR4 | IL6ST | MITF | NUTM1 | RANBP17 | SS18L1 | WIF1 |
| AXL | CDKN2A | EPHA3 | FH | IL7R | MLF1 | PALB2 | RAP1GDS1 | STAT3 | WISP3 |
| BAP1 | CDX2 | EPHA5 | FHIT | IRF4 | MLH1 | PAX3 | RARA | STAT4 | WRN |
| BARD1 | CHEK1 | EPHB1 | FIP1L1 | ITK | MLL1 | PAX5 | RB1 | STAT5B | WT1 |
| BCL10 | CHEK2 | EPH51 | FLCN | JAK1 | MLL10 | PAX7 | RBM15 | STIL | WWTR1 |
| BCL11A | CHIC2 | ERBB2 (HER2/NEU) | FLI1 | JAK2 | MLL13 | PBRM1 | REL | STK11 | XPA |
| BCL2L11 | CHN1 | ERBB3 (HER3) | FLT1 | JAK3 | MLL16 | PBX1 | RET | SUFU | XPC |
| BCL3 | CIC | ERBB4 (HER4) | FLT3 | JAZF1 | MN1 | PCM1 | RICTOR | SUZ12 | XPO1 |
| BCL6 | CIITA | ERC1 | FLT4 | KDMSA | MRE11 | PCSK7 | RMI2 | SYK | YWHAE |
| BCL7A | CLP1 | ERCC2 | FNBP1 | KDR (VEGFR2) | MSH2 | PDCC1 (PD1) | RNF43 | TAF15 | ZMYM2 |
| BCL9 | CLTC | ERCC3 | FOXA1 | KEAP1 | MSH6 | PDCC1LG2 (PDL2) | ROS1 | TCF12 | ZNF217 |
| BCR | CLTCL1 | ERCC4 | FOXO1 | KIAA1549 | MSI2 | PDGFB | RPL22 | TCF3 | ZNF331 |
| BIRC3 | CNBP | ERCC5 | FOXP1 | KIF5B | MTOR | PDGFRA | RPL5 | TCF7L2 | ZNF384 |
| BLM | CNTRL | ERG | FUBP1 | KIT | MYB | PDGFRB | RPN1 | TET1 | ZNF521 |
| BMPR1A | COPB1 | ESR1 | FUS | KLHL6 | MYC | PDK1 | RPTOR | TET2 | ZNF703 |
| BRAF | CREB1 | | | | | | | | |
| Whole Transcriptome Sequencing – Genes most commonly associated with cancer listed below. | | | | | | | | | |
| Fusions (RNA) | | | | | | | | Variant Transcripts (RNA) | |
| ABL | BRD3 | FGFR3 | INSR | MYB | NUMBL | PRKCA | RSPO3 | AR-V7 | |
| AKT3 | BRD4 | ERG | MAML2 | NOTCH1 | NUTM1 | PRKCB | TERT | | |
| ALK | EGFR | ESR1 | MAST1 | NOTCH2 | PDGFRA | RAF1 | TFE3 | | |
| ARHGAP26 | EWSR1 | ETV1 | MAST2 | NRG1 | PDGFRB | RELA | TFEB | EGFR vIII | |
| AXL | FGR | ETV4 | MET | NTRK1 | PIK3CA | RET | THADA | | |
| BCR | FGFR1 | ETV5 | MSMB | NTRK2 | PKN1 | ROS1 | TMPRSS2 | MET Exon 14 Skipping | |
| BRAF | FGFR2 | ETV6 | MUSK | NTRK3 | PPARG | RSPO2 | | | |

*Not available in New York State.

To order or learn more, visit www.Carismolecularintelligence.com.

US: 888.979.8669 | CustomerSupport@CarisLS.com – Intl: 00 41 21 533 53 00 | InternationalSupport@CarisLS.com

The Molecular Science Company

Caris Molecular Intelligence is a registered trademark of Caris MPI, Inc.

© 2019 Caris MPI, Inc. All rights reserved.

TN0221 v18 August 21, 2019

